
sanitas

SBF 14

DE	Diagnosewaage	
	Gebrauchsanweisung	2
EN	Diagnostic scale	
	Instruction for Use.....	10
FR	Pèse-personne impédancemètre	
	Mode d'emploi	17
IT	Bilancia diagnostica	
	Istruzioni per l'uso	25
RU	Диагностические весы	
	Инструкция по применению	33
NL	Diagnoseweegschaal	
	Gebruiksaanwijzing	41
PL	Waga diagnostyczna	
	Instrukcja obsługi	49

DE Deutsch

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Read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Contents

1. Included in delivery.....	10	6.1 Measuring weight only.....	13
2. Signs and symbols.....	10	6.2 Measuring your weight, carrying out diagnostics.....	13
3. Intended use.....	11	7. Evaluation of results.....	14
4. Warnings and safety notes.....	11	8. Incorrect measurement.....	15
5. Operation.....	13	9. Disposal.....	15
5.1 Inserting batteries.....	13	10. Warranty/service.....	16
5.2 Setting user data.....	13		
6. Operation.....	13		

1. Included in delivery


Check that the device packaging has not been tampered with and make sure that all contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

- SBF 14 diagnostic bathroom scale
- 2x 3V Lithium-battery CR 2032
- These instructions for use

2. Signs and symbols

The following symbols are used on the device, in these instructions for use, on the packaging and on the type plate for the device:

	WARNING Warning notice indicating a risk of injury or damage to health.
	IMPORTANT Safety note indicating possible damage to the device/accessory.
	Note Note on important information
	Observe the instructions for use
	Disposal in accordance with the Waste Electrical and Electronic Equipment EC Directive – WEEE
	Dispose of packaging in an environmentally friendly manner
	This product satisfies the requirements of the applicable European and national directives.
	The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.
	Manufacturer
	Certification symbol for products that are exported to the Russian Federation and members of the CIS.

 	Danger of tipping: position yourself in the centre of the weighing surface.	 	Danger of slipping: do not step on the scale with wet feet.
	Place the scale on an even surface. No carpet.		Do not overload the scale 180 kg / 396 lb / 28 st.

3. Intended use

The device is only intended for weighing humans and for recording your personal fitness data. The device is only intended for private use, and not for medical or commercial purposes.

4. Warnings and safety notes

WARNING

- **The scale may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.**
- Do not use during pregnancy.

- Important: do not step onto the scale with wet feet or if the surface of the scale is damp – danger of slipping!

- Do not step onto the outer edge of the scale on just one side: danger of tipping!

- Keep packaging material away from children (risk of suffocation).

Notes on handling batteries

- **CAUTION! RISK OF EXPLOSION!**

Never recharge, short circuit, mechanically deform, break up or open batteries. As a result, they may overheat, there is the risk of fire or they may split. Never use incorrect or different types of batteries. Never throw batteries into fire or water, or leave batteries on hot surfaces. The batteries may explode. Avoid

extreme conditions and temperatures above 50°C that could affect batteries, e.g. from radiators, fires, ovens. Otherwise there is a risk of explosion or leakage of flammable gas or liquid.

- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
-  **Choking hazard!** Small children may swallow and choke on batteries. Therefore, store batteries out of the reach of small children!
- Swallowing batteries can be extremely dangerous.
-  **Risk of chemical burns.**
 - If you think batteries have been swallowed or have otherwise entered the body, seek medical attention immediately.
 - If the battery compartment does not close securely, stop using the product and keep it away from children.
- Observe the plus (+) and minus (-) polarity signs.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- If the device is not to be used for a relatively long period, take the batteries out of the battery compartment.
- Do not use rechargeable batteries!
- Do not disassemble, open or crush the batteries.

General notes

- Please note that measuring tolerances are possible for technical reasons, as this scale is not calibrated for use in a professional medical context.
- The scale's maximum capacity is 180 kg (396 lb/28 st). The results for the weight measurement and bone mass calculation are displayed in 100-g increments (0.2 lb, 0.01 st).
- Measurements of body fat, body water and muscle content are displayed in increments of 0.1%.
- When supplied to the customer, the scale is set to weigh and measure in "kg" and "cm". To the rear of the scales, there is a switch where you can select "pound" and "stone" (lb, st).
- Place the scale on an even, hard surface; a hard surface is crucial for achieving accurate measurements.
- Repairs may only be carried out by Customer Services or authorised retailers. Before submitting a complaint, please check the batteries first and replace them if necessary.

Storage and maintenance

The accuracy of the measurements and service life of the device depend on its careful handling:

IMPORTANT

- The device should be cleaned from time to time. Do not use any abrasive cleaning products and never submerge the device in water.
- Make sure that no liquids come into contact with the scale. Never submerge the scale in water. Never rinse it in running water.
- Do not place any objects on the scale when it is not in use.
- Do not press the button violently or with pointed objects.
- Do not expose the scale to high temperatures or strong electromagnetic fields.
- Protect the device from knocks, damp, dust, chemicals, marked temperature fluctuations and nearby sources of heat (ovens, heaters).

General tips

- As far as possible, weigh yourself at the same time every day (ideally in the morning), after having been to the toilet, on an empty stomach and without clothing, in order to obtain results which can be compared.
- Important point regarding the measurement: the calculation of body fat may only be made when bare-foot; the soles of the feet may be lightly dampened where appropriate. Results may be unsatisfactory if the soles of the feet are completely dry or have a large amount of hard skin since conductivity will be impaired.
- Stand up straight and still during the measurement.
- Wait a few hours after any physical exertion to which your body is unaccustomed.
- Wait approx. 15 minutes after getting up so that the water stored in the body can disperse.
- Remember that only the long-term trend is important. Short-term changes in weight within a few days are normally caused by loss of fluids; however, body water plays an important role in terms of well-being.

Limitations

When calculating body fat and the other values, unusual and implausible results may occur in the case of:

- Children under approx. 10 years of age
- Performance athletes and body-builders
- Persons with a fever, persons undergoing dialysis, persons with oedema symptoms and persons suffering from osteoporosis
- Persons taking cardiovascular medicine (affecting the heart and vascular system)
- Persons taking vasodilative or vasoconstrictive medication
- Persons with significant anatomical abnormalities in their legs compared with their overall body height (leg length significantly shorter or longer).

5. Operation

5.1 Inserting batteries

If present, pull the battery insulating strip off the battery compartment cover or remove the battery’s protective film and insert the battery according to the polarity. If the scale fails to operate, remove the battery completely and insert it again. Your scale is equipped with a "replace battery" indicator. When the scale is operated with a battery which is too weak, "L0" appears in the display field and the scale automatically switches off. In this case the battery must be replaced (2 x 3V lithium battery CR2032).

5.2 Setting user data

"SET"

User data	Set values						
Memory position P	1 to 10						
Height	100 to 220 cm (3'-03" to 7'-03")						
Age	10-100						
Sex	male 🧑, female 🧒						
Degree of activity R	<table><tr><th>Physical activity</th></tr><tr><td>1 None.</td></tr><tr><td>2 Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).</td></tr><tr><td>3 Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.</td></tr><tr><td>4 High: physical exertion, at least 4 to 6 times a week for 30 minutes each.</td></tr><tr><td>5 Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.</td></tr></table>	Physical activity	1 None.	2 Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).	3 Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.	4 High: physical exertion, at least 4 to 6 times a week for 30 minutes each.	5 Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.
Physical activity							
1 None.							
2 Low: little and light physical exertion (e.g. going for walks, light gardening, simple gymnastics).							
3 Medium: physical exertion, at least 2 to 4 times a week for 30 minutes each.							
4 High: physical exertion, at least 4 to 6 times a week for 30 minutes each.							
5 Very high: intensive physical exertion, intensive training or hard physical work for at least 1 hour every day.							

Use the  and  buttons to select and confirm with "SET".

6. Operation

6.1 Measuring weight only

- Stand on the scale without moving about and distribute your weight evenly on both legs. The scale immediately begins to measure your weight. Your measured weight is displayed shortly after this (Fig. 1).
- The scale switches off a few seconds after you step down from it.

Fig. 1

6.2 Measuring your weight, carrying out diagnostics

- Switch on the scale (tap the standing surface).
- Select the memory preset where your basic personal data are stored by pressing the  or  button repeatedly. The data are displayed consecutively until "0.0" appears.
- Step onto the scale barefoot and make sure you're standing on both electrodes. First your body weight is determined and displayed (Fig. 1).
- Remain on the scale and avoid moving; your body fat is being analysed. This may take a few seconds.

The following data are displayed:

- Weight in kg (lb, st) (Fig. 1)
- Body fat percentage **FAT** in % (Fig. 2)
- Water percentage **BW** in % (Fig. 3)
- Muscle percentage **MUS** in % (Fig. 4)
- Bone mass in **BON** in kg (Fig. 5)
- Now all measured values are displayed consecutively and the scale switches off.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5

7. Evaluation of results

Body fat percentage

The following body fat percentages are guide values (contact your physician for further information).

Man

Age	low	normal	high	very high
10-14	<11	11-16	16,1-21	>21,1
15-19	<12	12-17	17,1-22	>22,1
20-29	<13	13-18	18,1-23	>23,1
30-39	<14	14-19	19,1-24	>24,1
40-49	<15	15-20	20,1-25	>25,1
50-59	<16	16-21	21,1-26	>26,1
60-69	<17	17-22	22,1-27	>27,1
70-100	<18	18-23	23,1-28	>28,1

Woman

Age	low	normal	high	very high
10-14	<16	16-21	21,1-26	>26,1
15-19	<17	17-22	22,1-27	>27,1
20-29	<18	18-23	23,1-28	>28,1
30-39	<19	19-24	24,1-29	>29,1
40-49	<20	20-25	25,1-30	>30,1
50-59	<21	21-26	26,1-31	>31,1
60-69	<22	22-27	27,1-32	>32,1
70-100	<23	23-28	28,1-33	>33,1

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values stated. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water percentage

The body water percentage is normally within the following ranges:

Man

Age	poor	good	very good
10-100	<50	50-65	>65

Woman

Age	low	normal	high
10-100	<65	45-60	>60

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage % is normally within the following ranges:

Man

Age	low	normal	high
10-14	<44	44-57	>57
15-19	<43	43-56	>56
20-29	<42	42-54	>54
30-39	<41	41-52	>52
40-49	<40	40-50	>50
50-59	<39	39-48	>48
60-69	<38	38-47	>47
70-100	<37	37-46	>46

Woman

Age	low	normal	high
10-14	<36	36-43	>43
15-19	<35	35-41	>41
20-29	<34	34-39	>39
30-39	<33	33-38	>38
40-49	<31	31-36	>36
50-59	<29	29-34	>34
60-69	<28	28-33	>33
70-100	<27	27-32	>32

Bone mass

Like the rest of our body, our bones are subject to natural growth, shrinking and ageing processes. The bone mass increases rapidly during childhood and reaches its peak between the ages of 30 and 40. As we age, our bone mass then begins to reduce. You can combat this reduction to an extent with the help of a healthy diet (particularly calcium and vitamin D) and regular physical exercise. The stability of the skeleton can be further increased by means of targeted muscle build-up. Please note that this scale does not identify the calcium content of the bones; instead, it calculates the weight of all the components that make up the bones (organic matter, inorganic matter and water). It is very difficult to influence the bone mass, although it does fluctuate within the scope of influencing factors (weight, height, age, gender). There are no recognised guidelines or recommendations.

IMPORTANT:

Please do not confuse bone mass with bone density.

The bone density can only be determined by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions on changes to bones and bone hardness (e.g. osteoporosis) using this scale.

8. Incorrect measurement

If the scale detects an error during weighing, “FFFF” or “Err” appears in the display.

If you step onto the scale before “0.0” appears in the display, the scale will not operate properly.

Possible causes of errors:	Remedy:
– The maximum load-bearing capacity of 180 kg (396 lbs) was exceeded.	– Only weigh the maximum permissible weight.
– The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	– Repeat weighing barefoot. – Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
– Your body fat lies outside the measurable range (less than 5% or greater than 55%).	– Repeat weighing barefoot. – Slightly moisten the soles of your feet if necessary.
– Scale is on carpeting.	– Place the scale on a level and solid surface.

9. Disposal

For environmental reasons, do not dispose of the device in the household waste at the end of its useful

life. Dispose of the unit at a suitable local collection or recycling point. Dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

Dispose of packaging in an environmentally friendly manner

The empty, completely flat batteries must be disposed of through specially designated collection boxes, recycling points or electronics retailers. You are legally required to dispose of the batteries.

The codes below are printed on batteries containing harmful substances:

Pb = Battery contains lead,

Cd = Battery contains cadmium,

Hg = Battery contains mercury.

10. Warranty/service

Hans Dinslage GmbH, Riedlinger Straße 28, 89524 Uttenweiler, Germany (hereinafter referred to as "HaDi") provides a warranty for this product, subject to the requirements below and to the extent described as follows.

The warranty conditions below shall not affect the seller's statutory warranty obligations which ensue from the sales agreement with the buyer.

The warranty shall apply without prejudice to any mandatory statutory provisions on liability.

HaDi guarantees the perfect functionality and completeness of this product.

The worldwide warranty period is 3 years, commencing from the purchase of the new, unused product from the seller.

The warranty only applies to products purchased by the buyer as a consumer and used exclusively for personal purposes in the context of domestic use.

German law shall apply.

During the warranty period, should this product prove to be incomplete or defective in functionality in accordance with the following provisions, HaDi shall carry out a repair or a replacement delivery free of charge, in accordance with these warranty conditions.

If the buyer wishes to make a warranty claim, they should approach their local retailer in the first instance: see the attached "International Service" list of service addresses.

The buyer will then receive further information about the processing of the warranty claim, e.g. where they can send the product and what documentation is required.

A warranty claim shall only be considered if the buyer can provide HaDi, or an authorised HaDi partner, with

- a copy of the invoice/purchase receipt, and
- the original product.

The following are explicitly excluded from this warranty:

- deterioration due to normal use or consumption of the product;
- accessories supplied with this product which are worn out or used up through proper use (e.g. batteries, rechargeable batteries, cuffs, seals, electrodes, light sources, attachments and nebuliser accessories);
- products that are used, cleaned, stored or maintained improperly and/or contrary to the provisions of the instructions for use, as well as products that have been opened, repaired or modified by the buyer or by a service centre not authorised by HaDi;
- damage that arises during transport between manufacturer and customer, or between service centre and customer;
- products purchased as seconds or as used goods;
- consequential damage arising from a fault in this product (however, in this case, claims may exist arising from product liability or other compulsory statutory liability provisions).

Repairs or an exchange in full do not extend the warranty period under any circumstances.

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