

sanitas

SFT 65

**SFT 65****Multifunktions-Thermometer**

6-in-1:

- Sekundenschnelle Messung im Ohr
- Sekundenschnelle Messung an der Stirn
- Oberflächentemperatur von Gegenständen und Flüssigkeiten
- Optischer Fieberalarm
- Datum und Uhrzeit
- 10 Speicherplätze

Geeignet für die ganze Familie

Einfache Reinigung

Großes, einfach ablesbares Display

Anzeige von normaler oder erhöhter Temperatur durch farbige LED

- grün = alles in Ordnung
- rot = Fieberalarm ($\geq 38^{\circ}\text{C}$ / $\geq 100,4^{\circ}\text{F}$)

Messgenauigkeit:

- Ohr/Stirn: $\pm 0,2^{\circ}\text{C}$ bei $35 - 42^{\circ}\text{C}$
($\pm 0,4^{\circ}\text{F}$ bei $95 - 107,6^{\circ}\text{F}$)
- Im restlichen Messbereich: $\pm 0,3^{\circ}\text{C}$ ($\pm 0,5^{\circ}\text{F}$)
- Objekt: $\pm 1,5^{\circ}\text{C}$ ($\pm 2,7^{\circ}\text{F}$) bei $<30^{\circ}\text{C}$;
 $\pm 5\%$ bei $>30^{\circ}\text{C}$

Anzeige von $^{\circ}\text{C}$ und $^{\circ}\text{F}$

Moderne, handliche Optik

Abschaltautomatik

Mit Batteriewechselanzeige

Inklusive 2 x 1,5V Batterie AAA

Produktmaße (LxBxH): 138 x 35 x 33 mm

Produktgewicht: ca. 90 g

Medizinprodukt

2 Jahre Garantie

VE: 4 / Transportkarton: 8 x 4

EAN-Nr.: 42 11125 79515 3

Art.-Nr.: 795.15

SFT 65**Multi-function thermometer**

6-in-1:

- Measurement within seconds in the ear
- Forehead temperature measurement in seconds
- Surface temperature of objects and liquids
- Optical high temperature alarm
- Date and time
- 10 memory spaces

Suitable for the whole family

Easy to clean

Large, easy to read display

Indication of normal or heightened temperature via coloured LED

- green = everything o.k.
- red = fever alarm ($\geq 38^{\circ}\text{C}$ / $\geq 100,4^{\circ}\text{F}$)

Measurement accuracy:

- Ear/Forehead: $\pm 0,2^{\circ}\text{C}$ with $35 - 42^{\circ}\text{C}$
($\pm 0,4^{\circ}\text{F}$ with $95 - 107,6^{\circ}\text{F}$)
- In remaining measurement range: $\pm 0,3^{\circ}\text{C}$ ($\pm 0,5^{\circ}\text{F}$)
- Object: $\pm 1,5^{\circ}\text{C}$ ($\pm 2,7^{\circ}\text{F}$) with $<30^{\circ}\text{C}$;
 $\pm 5\%$ with $>30^{\circ}\text{C}$

Indication of $^{\circ}\text{C}$ and $^{\circ}\text{F}$

Modern, handy optic

Automatic switch off

Mit Batteriewechselanzeige

Including 2 x 1,5V Battery AAA

Product measurements (LxWxH): 138 x 35 x 33 mm

Product weight: approx. 90 g

Medical device

2 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 42 11125 79515 3

Item no.: 795.15