

# sanitas

Fieber-  
thermometer



SFT 53

## SFT 53 Ohrthermometer

Hygienisch und sicher messen

Mit auswechselbarer Schutzkappe

Infrarot-Messtechnik

Sekundenschnelle Messung

Fieberalarm ab 38°C

Warnton ab 37,5°C

Anzeige von °C und °F

Objekttemperaturmessung

Anzeige der Raumtemperatur

10 Speicherplätze

Signalton bei Messende

Lachendes Smiley = alles in Ordnung

Trauriges Smiley =  $\geq 38^{\circ}\text{C}$  Fieberalarm

Abschaltautomatik

Mit Batteriewechselanzeige

Mit Aufbewahrungbox

Inkl. 10 Schutzkappen

Inkl. 2 x 1,5V Batterien AAA

Produktmaße: 155 x 38 x 47 mm

Produktgewicht: ca. 63 g (ohne Batterien)

2 Jahre Garantie

### Nachkaufset SFT 53 Schutzkappen

20 Schutzkappen (2 x 10)

VE: 10 / Transportkarton: 10 x 10

EAN-Nr.: 4211125 79526 9

Art.-Nr.: 795.26

VE: 4 / Transportkarton: 12 x 4

EAN-Nr.: 4211125 79525 2

Art.-Nr.: 795.25

## SFT 53 Digital thermometer

Hygienic and safe measurement

With interchangeable protective cap

Infrared measurement technology

Measurement in seconds

Temperature alarm at 38°C

Warning sound at 37,5°C

Indication of °C and °F

Object temperature measurement

Room temperature display

10 memory spaces

Acoustic signal when measurement is complete

Laughing Smiley = everything is OK

Sad Smiley =  $\geq 38^{\circ}\text{C}$  temperature alarm

Automatic switch off

With low battery indicator

With storage box

Incl. 10 protective caps

Incl. 2 x 1.5V batteries AAA

Product measurements: 155 x 38 x 47 mm

Product weight: approx. 63 g (without batteries)

2 year guarantee

### Replacement Set SFT 53 Protective Caps

20 Protective caps (2 x 10)

Sales unit: 10 / Shipping carton: 10 x 10

EAN no.: 4211125 79526 9

Item no.: 795.26

Sales unit: 4 / Shipping carton: 12 x 4

EAN no.: 4211125 79525 2

Item no.: 795.25