



SANITAS

SBF 14



ETM TESTMAGAZIN · URTEIL	
Sanitas SBF 14	
GUT	89,21 %
Glas-Diagnosewaage im Onlinetest www.etm-testmagazin.de	04/2016

SBF 14 Glas-Diagnosewaage

Quadratische Trittpläche aus Sicherheitsglas
Moderne Sensortasten
Gut lesbare LCD-Anzeige: Ziffernhöhe 26 mm
Elektroden aus gebürstetem Edelstahl
Ermittlung von Körpergewicht, -fett, -wasser, Muskelanteil und Knochenmasse
Einteilung: - 0,1% für Körperfett, -wasser und Muskelanteil - 100 g für Körpergewicht und Knochenmasse
5 Aktivitätsgrade
10 Benutzerspeicher
Umstellung von kg, lb, st
Tragkraft: 180 kg
Überlastungsanzeige
Abschaltautomatik
Einschalttechnik: Step-On
Inklusive 2 x 3V CR2032 Batterie
Produktmaße (LxBxH): 30 x 30 x 2,3 cm
Produktgewicht: ca. 1550 g
3 Jahre Garantie
VE: 4 / Transportkarton: -
EAN-Nr.: 4211125 73591 3
Art.Nr.: 735.91

SBF 14 Glass diagnostic scale

Square platform of safety glass
Modern sensor buttons
Easy to read LCD-display: digit size 26 mm
Electrodes of brushed stainless steel
Display of body weight, body fat, body water, muscle ratio and bone mass
Graduation: - 0.1% for body fat, body water and muscle ratio - 100 g for body weight and bone mass
5 activity levels
10 user memories
Change over to kg, lb, st
Capacity: 180 kg
Overload indication
Automatic switch-off
Switch-on technology: Step-On
2 x 3V CR2032 battery included
Product measurements (LxWxH): 30 x 30 x 2.3 cm
Product weight: approx. 1550 g
3 year guarantee
Sales unit: 4 / Shipping unit: -
EAN no.: 4211125 73591 3
Item no.: 735.91