

SANITAS

SEM 43



3 in 1:

Muskeltraining und
-regeneration durch **EMS**
(Elektrische Muskelstimulation)



Einfache Schmerzlinderung
durch **TENS** (Transkutane
Elektrische Nervenstimulation)



Entspannung und
Massage



SEM 43 Digital EMS/TENS

Funktionen:

- Elektrische Muskelstimulation (EMS) - für Muskeltraining und -regeneration
- Transkutane Elektrische Nervenstimulation (TENS) - zur medikamenten- und nebenwirkungsfreien Schmerzlinderung
- Massage - für Entspannung und Wohlbefinden

2 getrennt regelbare Kanäle mit 4 selbstklebenden Elektroden

44 vorprogrammierte Anwendungen (EMS/TENS/Massage)

6 individualisierbare Programme (EMS/TENS)

"Doctor's-Function": zur idealen Abstimmung auf Ihr persönliches Therapieprogramm

Timer-Funktion (5-100 Minuten einstellbar)

Tastensperre

Hintergrundbeleuchtetes, leicht lesbares Display

Ergonomisches Design, mit Gürtelclip

Sicherheitsabschaltung

Ausgangsstrom: max. 200 mA p-p (500 Ohm)

Ausgangsfrequenz: 1-150 Hz

Pulsbreite: 50-450 µs pro Puls

Impulsform: biphasische Rechteckimpulse

Mit ausführlicher Gebrauchsanweisung (inkl. Anwendungsmöglichkeiten)

Medizinprodukt

Inklusive 3 x AAA Batterien

Produktmaße: 135 x 66 x 29 mm

Produktgewicht: ca. 107 g

2 Jahre Garantie

Versandereinheit: 4

Transportkarton: 6 x 4

EAN-Nr.: 4211125 66130 4

Art.-Nr.: 661.30

SEM 43 Digital EMS/TENS

Functions:

- Electronic muscle stimulation (EMS) - for muscle training and regeneration
- Transcutaneous electrical nerve stimulation (TENS) - medicine- and side effect-free pain-relief
- Massage - for relaxation and well-being

2 separately adjustable channels with 4 self-adhesive electrodes

44 pre-programmed applications (EMS/TENS/Massage)

6 individualizable programs (EMS/TENS)

"Doctor's-Function": ensuring ideal coordination to your personal therapy programme

Timer function (adjustable between 5-100 minutes)

Key lock

Backlight illuminated, easy readable display

Ergonomic design, with belt clip

Safety switch-off

Output current: max. 200 mA p-p (500 Ohm load)

Output frequency: 1-150 Hz

Pulse width: 50-450 µs per pulse

Pulse form: dual phase square pulses

With detailed instruction manual (inclusive application possibilities)

Medical product

3 x AAA batteries included

Product measurements: 135 x 66 x 29 mm

Product weight: approx. 107 g

2 year guarantee

Sales unit: 4

Shipping carton: 6 x 4

EAN no.: 4211125 66130 4

Item no.: 661.30