

# SANITAS



SFT 08

## SFT 08

### Thermometer

Digitales Fieberthermometer

Messung in ca. 30 Sekunden

Flexible Mess-Spitze

Ohne Glas, quecksilberfrei

1 Speicherplatz

Messgenauigkeit  $\pm 0,1^{\circ}\text{C}$   
(zwischen  $35,5$  und  $42^{\circ}\text{C}$ )

Messgenauigkeit  $\pm 0,2^{\circ}\text{C}$   
(zwischen  $32 - 35,4^{\circ}\text{C}$  und  $42,1 - 42,9^{\circ}$ )

Signalton bei Messende

Warnton, wenn Messergebnis über  $37,8^{\circ}\text{C}$   
(Fieberalarm)

Extra großes Display

Wasserdicht

Desinfektionsfähig

Mit Schutzkappe

Batteriewechselanzeige

Abschaltautomatik

Medizinprodukt

Lieferung im praktischen Verkaufsdisplay

Anzeige von  $^{\circ}\text{C}$  und  $^{\circ}\text{F}$

Inklusive 1x Batterie LR41

Produktmaße (LxBxH): 142 x 23 x 13 mm

Produktgewicht: ca. 14 g (mit Batterien)

2 Jahre Garantie

VE: 20 / Transportkarton: 10 x 20

EAN-Nr.: 4211125 79108 7

Art.-Nr.: 791.08

## SFT 08

### Thermometer

Digitales thermometer

Measurement in approx. 30 seconds

Flexible measuring tip

Without glass, mercury free

1 memory

Measurement accuracy  $\pm 0.1^{\circ}\text{C}$   
(between  $35.5$  and  $42^{\circ}\text{C}$ )

Measurement accuracy  $\pm 0.2^{\circ}\text{C}$   
(between  $32 - 35.4^{\circ}\text{C}$  and  $42.1 - 42.9^{\circ}\text{C}$ )

Audible signal when measuring stops

Audible warning signal if temperature exceeds  $37,8^{\circ}\text{C}$   
(fever alarm)

Extra large display

Waterproof

Can be disinfected

With protection cap

Battery changing indicator

Automatic switch off

Medical product

Delivered in practical sales display

Indication of  $^{\circ}\text{C}$  and  $^{\circ}\text{F}$

Including 1x Battery LR41

Product measurements (LxWxH): 142 x 23 x 13 mm

Product weight: approx. 14 g (with batteries)

2 year guarantee

Sales unit: 20 / Shipping carton: 10 x 20

EAN no.: 4211125 79108 7

Item no.: 791.08