


SANITAS



SHK 32



## SHK 32

### Rücken-/Nacken Heizkissen

Superflauschige Oberfläche, atmungsaktiv

Anpassungsfähig und hautsympatisch

Individuell einstellbar durch Bauchgurt und Klettverschluss am Hals

3 beleuchtete Temperaturstufen

Abschaltautomatik (nach ca. 90 Minuten)

Maschinenwäsche bei 30°C

Abnehmbarer Schalter

Sicherheitssystem

Elektronische Temperatur-Regelung

Schnellheizung

Öko-Tex 100

Die bei diesem Gerät eingesetzten Textilien erfüllen die hohen humanökologischen Anforderungen des Öko-Tex Standards 100, wie durch das Forschungsinstitut Hohenstein nachgewiesen.

100 Watt

Produktmaße: 54 x 38 cm

Produktgewicht: ca. 390 g

2 Jahre Garantie

VE: 4 / Transportkarton: -

EAN-Nr.: 4211125 21516 3

Art.-Nr.: 215.16

## SHK 32

### Back-/Neck heating pad

Super soft surface, breathable

Cosy and skinfriendly

Individually adjustable using abdominal belt and hook-and-loop fastener at neck

3 illuminated temperature settings

Automatic switch off (after approx. 90 min.)

Machine washable at 30°C

Detachable control

Safety system

Electronic temperature control

Fast heat up

Oeko-Tex 100

The textiles used for this device meet the stringent human ecological requirements of Oeko-Tex Standard 100, as verified by Hohenstein Research Institute.

100 watt

Product measurements: 54 x 38 cm

Product weight: approx. 390 g

2 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 4211125 21516 3

Item no.: 215.16