

# sanitas



## INHALATOR | NEBULISER

SIH 09



Verneblung mit  
Kompressor-Druck-  
luft-Technologie



Umfangreiches  
Zubehör



Zur Behandlung der oberen und unteren Atemwege

Verneblung von flüssigen Medikamenten mittels  
Kompressor-Drucklufttechnologie

Zur Anwendung bei Erkältung, Asthma,  
Atemwegserkrankung ...

Hohe Verneblungsleistung (ca. 0,3 ml/min)  
- kurze Inhalationszeit

Zubehör: Mundstück, Erwachsenen und Kindermaske inkl.  
Halteband, Vernebler, 5 Filter, Schlauch

Desinfektionsfähig

0,35 - 0,8 bar Arbeitsdruck

230 V~ / 50 Hz / 230 VA

Produktmaße (LxBxH): 133 x 106 x 90 mm

Produktgewicht: ca. 0,85 ±0,2 kg

Medizinprodukt

3 Jahre Garantie

For the treatment of the upper and lower respiratory tract

Nebulisation of liquid medicine using compressed  
air technology

For use in colds, asthma, respiratory disease ...

High nebulisation capacity (approx. 0.3 ml/min)  
- short inhalation time

Accessories: mouthpiece, adult mask and children's mask  
including strap, atomizer, 5 filters, hose

Can be disinfected

0.35 - 0.8 bar operating pressure

230 V~ / 50 Hz / 230 VA

Product measurements (LxWxH): 133 x 106 x 90 mm

Product weight: approx. 0.85 ±0.2 kg

Medical product

3 year guarantee

