


sanitas

Diagnose-
waage



SBF 14

SBF 14 Diagnosewaage

Quadratische Trittfläche aus Sicherheitsglas

Moderne Sensortasten

Gut lesbare LCD-Anzeige: Ziffernhöhe 26 mm

Elektroden aus gebürstetem Edelstahl

Ermittlung von Körpergewicht, -fett, -wasser, Muskelanteil und Knochenmasse

Einteilung:

- 0,1% für Körperfett, -wasser und Muskelanteil
- 100 g für Körpergewicht und Knochenmasse

5 Aktivitätsgrade

10 Benutzerspeicher

Umstellung von kg, lb, st

Tragkraft: 180 kg

Überlastungsanzeige

Abschaltautomatik

Einschalttechnik: Step-On

Inklusive 2 x 3V CR2032 Batterie

Produktmaße (LxBxH): 30 x 30 x 2,3 cm

Produktgewicht: ca. 1550 g

3 Jahre Garantie

VE: 4 / Transportkarton: -

EAN-Nr.: 4211125 73591 3

Art.Nr.: 735.91

SBF 14 Diagnostic scale

Square platform of safety glass

Modern sensor buttons

Easy to read LCD-display: digit size 26 mm

Electrodes of brushed stainless steel

Display of body weight, body fat, body water, muscle ratio and bone mass

Graduation:

- 0.1% for body fat, body water and muscle ratio
- 100 g for body weight and bone mass

5 activity levels

10 user memories

Change over to kg, lb, st

Capacity: 180 kg

Overload indication

Automatic switch-off

Switch-on technology: Step-On

2 x 3V CR2032 battery included

Product measurements (LxWxH): 30 x 30 x 2.3 cm

Product weight: approx. 1550 g

3 year guarantee

Sales unit: 4 / Shipping unit: -

EAN no.: 4211125 73591 3

Item no.: 735.91