


sanitas

Diagnose-
waage



SBG 21



plus +

Alles auf
einen Blick

SBG 21 Diagnosewaage

Grafisches Druckmotiv mit „Durchblick“

Großes LCD-Display mit Multi-Anzeige

Ziffernhöhe (Gewichtsanzeige): 25 mm

Körpergewicht, -fett, -wasser, Muskelanteil,
Knochenmasse, Kalorienbedarf (BMR/AMR)

Mit BMI-Berechnung

Interpretation des Körperfetts

Elektroden aus gebürstetem Edelstahl

5 Aktivitätsgrade

10 Speicherplätze

Überlastungsanzeige

180 kg Tragkraft

Einteilung:

- 0,1% für Körperfett, -wasser, Muskelanteil
- 100 g für Gewicht, Knochenmasse

Tap-On-Technik

Abschaltautomatik

Inkl. 2 x 3 V Batterie CR2032

Produktabmessungen (LxWxH): 30 x 31,5 x 2,5 cm

Produktgewicht: ca. 1600 g

3 Jahre Garantie

Versandereinheit: 4 / Transportkarton: -

EAN-Nr.: 42 11125 76435 7

Art.-Nr.: 764.35

SBG 21 Diagnostic bathroom scale

Printing with „see-through“-effect

Large LCD-display with multi display

Digit size (weight display): 25 mm

Body weight, body fat, body water, muscle ratio,
bone mass, calory consumption (BMR/AMR)

With BMI calculation

Interpretation of body fat

Electrodes of brushed stainless steel

5 activity levels

10 memories

Overload indication

180 kg capacity

Division:

- 0,1% body fat, body water, muscle ratio
- 100 g for body weight, bone mass

Tap-on-technology

Automatic switch off

2 x 3 V batteries CR2032 included

Product measurements (LxWxH): 30 x 31.5 x 2.5 cm

Product weight: approx. 1600 g

3 year guarantee

Sales unit: 4 / Shipping carton: -

EAN no.: 42 11125 76435 7

Item number: 764.35