

sanitas

INHALATOR | NEBULISER



Verneblung mit
Kompressor-Druck-
luft-Technologie



SIH 21



Zur Behandlung der oberen und unteren Atemwege

Verneblung von flüssigen Medikamenten mittels
Kompressor-Drucklufttechnologie

Zur Anwendung bei Erkältung, Asthma,
Atemwegserkrankung ...

Hohe Vernebelungsleistung (ca. 0,4 ml/min)
- kurze Inhalationszeit

Zubehör: Mundstück, Erwachsenenmaske, Kindermaske,
Druckluft-Schlauch, Filter

Mit praktischem Aufbewahrungsfach für die Zubehörteile
Desinfektionsfähig

0,5 - 1,45 bar Arbeitsdruck

230 V~ / 50 Hz / 230 VA

Produktmaße (LxBxH): 295 x 177 x 97 mm

Produktgewicht: ca. 2,1 kg

Medizinprodukt

2 Jahre Garantie

Typ: SIH 21/2

VE: 2

EAN-Nr.: 4211125 60114 0

Art.-Nr.: 601.14

For the treatment of the upper and lower respiratory tract

Nebulisation of liquid medicine using compressed
air technology

For use in colds, asthma, respiratory disease ...

High nebulisation capacity (approx. 0.4 ml/min)
- short inhalation time

Accessories: Mouthpiece, adult mask, child mask,
compressed air tube, filter

With practical storage compartment for the accessories
Can be disinfected

0.5 - 1.45 bar operating pressure

230 V~ / 50 Hz / 230 VA

Product measurements (LxWxH): 295 x 177 x 97 mm

Product weight: approx. 2.1 kg

Medical product

2 year guarantee

Type: SIH 21/2

Sales unit: 2

EAN no.: 4211125 60114 0

Item no.: 601.14