

sanitas

SFT 08



SFT 08 Thermometer

Digitales Fieberthermometer
Messung in ca. 30 Sekunden
Flexible Mess-Spitze
Ohne Glas, quecksilberfrei
1 Speicherplatz
Messgenauigkeit $\pm 0,1^{\circ}\text{C}$ (zwischen $35,5$ und 42°C)
Messgenauigkeit $\pm 0,2^{\circ}\text{C}$ (zwischen $32 - 35,4^{\circ}\text{C}$ und $42,1 - 42,9^{\circ}\text{C}$)
Signalton bei Messende
Warnton, wenn Messergebnis über $37,8^{\circ}\text{C}$ (Fiebertalarm)
Extra großes Display
Wasserdicht
Desinfektionsfähig
Mit Schutzkappe
Batteriewechselanzeige
Abschaltautomatik
Medizinprodukt
Lieferung im praktischen Verkaufsdisplay
Anzeige von $^{\circ}\text{C}$ und $^{\circ}\text{F}$
Inklusive 1x Batterie LR41
Produktmaße (LxBxH): $142 \times 23 \times 13 \text{ mm}$
Produktgewicht: ca. 14 g (mit Batterien)

2 Jahre Garantie

VE: 20 / Transportkarton: 10×20
EAN-Nr.: 4211125 79108 7
Art.-Nr.: 791.08

SFT 08 Thermometer

Digital thermometer
Measurement in approx. 30 seconds
Flexible measuring tip
Without glass, mercury free
1 memory
Measurement accuracy $\pm 0.1^{\circ}\text{C}$ (between 35.5 and 42°C)
Measurement accuracy $\pm 0.2^{\circ}\text{C}$ (between $32 - 35.4^{\circ}\text{C}$ and $42.1 - 42.9^{\circ}\text{C}$)
Audible signal when measuring stops
Audible warning signal if temperature exceeds $37,8^{\circ}\text{C}$ (fever alarm)
Extra large display
Waterproof
Can be disinfected
With protection cap
Battery changing indicator
Automatic switch off
Medical product
Delivered in practical sales display
Indication of $^{\circ}\text{C}$ and $^{\circ}\text{F}$
Including 1x Battery LR41
Product measurements (LxWxH): $142 \times 23 \times 13 \text{ mm}$
Product weight: approx. 14 g (with batteries)

2 year guarantee

Sales unit: 20 / Shipping carton: 10×20
EAN no.: 4211125 79108 7
Item no.: 791.08